Tell us whether you accept cookies

We use <u>cookies to collect information</u> about how you use GOV.UK. We use this information to make the website work as well as possible and improve government services.

Accept all cookies

Set cookie preferences





Coronavirus (COVID-19): what you need to do

Stay at home

- Only go outside for food, health reasons or work (but only if you cannot work from home)
- If you go out, stay 2 metres (6ft) away from other people at all times
- · Wash your hands as soon as you get home

Do not meet others, even friends or family.

You can spread the virus even if you don't have symptoms.

Hide message

Home > COVID-19: guidance for sampling and for diagnostic laboratories



Guidance

COVID-19: guidance on alternative swab types and transport media

Updated 16 April 2020

Difficulties with supply chain for key virological testing components such as swabs and transport medium have led to local variations in practice with respect to sampling which may also impact on PCR diagnostic test performance.

In addition, it is important that we make best use of our available testing capacity by utilising a single swab wherever practical. To mitigate the supply issues, recommendations are two-fold:

- 1.Reduce demand for swab usage by:
- · using a single swab for sampling throat then nose
- reduce level of testing for non-COVID related viruses (for example influenza and other respiratory viruses, as perRoyal College of Pathologists guidance

2.Use alternative swabs and transport medium in accordance with a locally validated laboratory strategy for preferred and secondary source of both swabs and transport medium. Recommendations for alternatives are provided below.

Alternative swab types

Multiple reports indicate that flocked or rayon or Dacron swabs are the preferred types for collecting respiratory samples, with decreased sensitivity with cotton swabs. Urethral, pernasal or mini-tip swabs may also be alternatives if standard sized swabs are not available. Dry swabs are not optimal for respiratory samples but can be used if no alternative is available. 'Charcoal' swabs (such as swabs transported in a charcoal-containing medium) are not suitable.

Swabs with polystyrene shafts should be used, as swabs with wooden shafts may contain PCR inhibitors and may also be a safety issue, especially for patients who are self-swabbing or when used to swab young children.

Alternative transport medium

Laboratories can use alternatives to Viral Transport Medium (VTM) subject to local validation, including Universal Transport Medium (UTM) or liquid Amies medium without charcoal (Amies medium with charcoal is not suitable). If necessary the following can be used as alternatives (in order of hierarchy): Phosphate Buffered Saline (PBS) with a protein source such as albumin or bovine serum, PBS alone, saline, or no transport medium (such as a dry swab).

Is this page useful?

Yes

No

Is there anything wrong with this page?

Coronavirus (COVID-19)

Coronavirus (COVID-19): what you need to do

Transition period

Transition period: check how to get ready

Services and information

Benefits

Births, deaths, marriages and care

Business and self-employed

Childcare and parenting

Citizenship and living in the UK

Crime, justice and the law

Disabled people

Driving and transport

Education and learning

Employing people

Environment and countryside

Housing and local services

Money and tax

Passports, travel and living abroad

Visas and immigration

Working, jobs and pensions

Departments and policy

How government works

Departments

Worldwide

Services

Guidance and regulation

News and communications

Research and statistics

Policy papers and consultations

Transparency and freedom of information releases

<u>Help Privacy Cookies Contact Accessibility statement Terms and conditions Rhestr o Wasanaethau Cymraeg</u>
Built by the <u>Government Digital Service</u>

OGL

All content is available under the Open Government Licence v3.0, except where otherwise stated



© Crown copyright